



don whitley
scientific
excellence in microbiology

Sample Preparation easyMIX Blender



New low-price blender

**easyMIX is a quality blender
at an attractive price**

easyMIX is simple to use – the programmable timer allows you to select blending cycles from ten seconds to three minutes or a continuous mode.

The blending cycle automatically starts when the door is closed. easyMIX has a removable door, removable paddles and a wipe-down keypad so is easy to keep clean and an optional waste drawer contains any spills.

easyMIX is a quality instrument at a competitive price – the chamber and blending paddles are manufactured from stainless steel.

Quiet in operation, the easyMIX is suitable for 80ml* and 400ml blending bags.

**Specific positioning required*

www.dwscientific.co.uk

Standard Features

Digital control
Programmable timer: 10 secs to 3 mins or continuous mode
Automatic start of cycle upon door closing
Wipe clean keypad
Easy-to-clean, fully opening door, removable paddles

Options

S01972 Waste drawer

Accessories

S00172 Blending bag rack to hold 10 bags (perspex)

S01744 Blender Bag Holder/Opener (holder only)

S01741 Rack to hold 10 blender bags (acrylic)

B00199 Bag closing clips (pack of 50)



Specification

Order Code: S02125
Power supply: 230V - 50/60Hz
Instrument Size : 600mm x 400mm x 275mm (w x d x h)
Instrument Weight: 18kg

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved.
© 2006 Don Whitley Scientific Limited.

Don Whitley Scientific Limited

14 Otley Road, Shipley, West Yorkshire, BD17 7SE, England.

Telephone: +44 (0)1274 595728 Fax: +44 (0)1274 531197

Website: www.dwscientific.co.uk Email: info@dwscientific.co.uk