

THEY THINK OF EVERYTHING!

DON WHITLEY SCIENTIFIC INTRODUCES AUTOMATIC PLATE LABELLING

When it comes to sample preparation, Don Whitley Scientific have got it covered - from start to finish. The latest addition to their range, the AES Labelmaster, sees the introduction of automatic Petri dish labelling, providing an accurate and reliable alternative to the tedious and time-consuming manual method.

The new system, which works equally well with full or empty dishes, optimises the preparation and follow-up of sample analysis and reduces mistakes due to human error. With an impressive capacity of 500 plates and a production rate of up to 1200 plates an hour, the speedy Labelmaster is ideal for large contract laboratories processing thousands of plates a week.

Software is included for use on a stand-alone basis, but the system is also compatible with the AES APS300 pourer stacker, the AES Labmaster, and LIMS.

Users may choose from the following parameters when setting up a label: sample or patient name, medium, bacteria, dilution, sample volume, batch number and bar code – turning a chore into a couple of minutes' work.

Editor: for further information, please contact Gill Iredale at Don Whitley Scientific, 14 Otley Road, Shipley, West Yorkshire, BD17 7SE.

Tel: 01274 595728. Fax: 01274 531197

email: info@dwscientific.co.uk