



don whitley
scientific
excellence in microbiology

Whitley Workstations Sleeve Kits



Supplied as ready-assembled units

Individual components can be ordered

Neoprene cuffs are available for anyone concerned about possible latex allergies

Whitley Workstation sleeve systems consist of components produced in different materials to offer the best possible combination of system longevity, comfort and freedom of movement for the operator, and an effective seal to the workstation.

Wrist seals are available in either latex or neoprene. Sleeves are manufactured from polyurethane. The systems are available either as ready-assembled kits or as individual components.

We continue to sell a range of latex gauntlets in four sizes.

Notes

Latex wrist seals can be trimmed to fit larger wrists. Trim the seals using sharp scissors or a scalpel, taking care not to tear the material or leave nicks from which a tear could start.

Neoprene is naturally fragile and care should be exercised during handling and use. Neoprene seals should not be trimmed.



Complete Sleeve Kits

Comprises: two wrist seals, two 'o' rings, two cuff rings and two sleeves

A02820	Pair of Sleeve Kits with latex wrist seals (one size fits all)
A02830	Pair of Sleeve Kits with Neoprene wrist seals (standard size)
A02840	Pair of Sleeve Kits with Neoprene wrist seals (large size)

Replacement Components for Latex Sleeve Kits

A02821	Pair of latex wrist seals (one size fits all)
A02715	White Polyurethane Sleeve (<i>single</i>)
A02716	Sleeve cuff ring (<i>single</i>)
A02672	Sleeve retaining 'O' ring (<i>single</i>)

Replacement Components for Neoprene Sleeve Kits

A02831	Pair of Neoprene wrist seals (standard size)
A02841	Pair of Neoprene wrist seals (large size)
A02715	White Polyurethane Sleeve (<i>single</i>)
A02716	Sleeve cuff ring (<i>single</i>)
A02672	Sleeve securing 'O' ring

Gauntlets - not illustrated

A02433	Pair of latex gauntlets (size 5)
A02645	Pair of latex gauntlets (size 6)
A02647	Pair of latex gauntlets (size 8)
A02648	Pair of latex gauntlets (size 9)

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved. © 2007 Don Whitley Scientific Limited.

Don Whitley Scientific Limited

14 Otley Road, Shipley, West Yorkshire, BD17 7SE, England.

Telephone: +44 (0)1274 595728 Fax: +44 (0)1274 531197

Website: www.dwsscientific.co.uk Email: info@dwsscientific.co.uk