



don whitley  
scientific  
excellence in microbiology

## Whitley Workstations Sleeve Kits



**Supplied as ready-assembled units**

**Individual components can be ordered**

**Neoprene cuffs are available for anyone concerned about possible latex allergies**

Whitley Workstation sleeve systems consist of components produced in different materials to offer the best possible combination of system longevity, comfort and freedom of movement for the operator, and an effective seal to the workstation.

Wrist seals are available in either latex or neoprene. Sleeves are manufactured from polyurethane. The systems are available either as ready-assembled kits or as individual components.

We continue to sell a range of latex gauntlets in three sizes.

## Notes

Latex wrist seals can be trimmed to fit larger wrists. Trim the seals using sharp scissors or a scalpel, taking care not to tear the material or leave nicks from which a tear could start.

Neoprene is naturally fragile and care should be exercised during handling and use. Neoprene seals should not be trimmed.



## Complete Sleeve Kits

two x wrist seals, two x 'o' rings, two x cuff rings and two x sleeves

**A02820** Pair of Sleeve Kits with latex wrist seals (one size fits all)

**A02830** Pair of Sleeve Kits with Neoprene wrist seals (standard size)

**A02840** Pair of Sleeve Kits with Neoprene wrist seals (large size)

## Replacement Components for Latex Sleeve Kits

**A02821** Pair of latex wrist seals (one size fits all)

**A02715** White Polyurethane Sleeve (*single*)

**A02716** Sleeve cuff ring (*single*)

## Replacement Components for Neoprene Sleeve Kits

**A02831** Pair of Neoprene wrist seals (standard size)

**A02841** Pair of Neoprene wrist seals (large size)

**A02715** White Polyurethane Sleeve (*single*)

**A02716** Sleeve cuff ring (*single*)

## Gauntlets - not illustrated

**A00121** Pair of latex gauntlets (size 6)

**A00118** Pair of latex gauntlets (size 8)

**A00123** Pair of latex gauntlets (size 10)

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved. © 2007 Don Whitley Scientific Limited.

## Don Whitley Scientific Limited

14 Otley Road, Shipley, West Yorkshire, BD17 7SE, England.

Telephone: +44 (0)1274 595728 Fax: +44 (0)1274 531197

Website: [www.dwscientific.co.uk](http://www.dwscientific.co.uk) Email: [info@dwscientific.co.uk](mailto:info@dwscientific.co.uk)