



don whitley
scientific

PRESS RELEASE

Don Whitley Scientific Limited
14 Otley Road, Shipley
West Yorkshire, BD17 7SE, UK
Tel: +44 (0) 1274 595728
Fax: +44 (0) 1274 531197
E-mail: info@dwscientific.co.uk
Web: www.dwscientific.co.uk

A85 – MAXIMUM FLEXIBILITY TO SUIT YOUR ANAEROBIC WORKLOAD

Don Whitley Scientific has unveiled the new Whitley A85 Anaerobic workstation to offer increased workload flexibility for microbiology and cell culture applications. The innovative A85 design features patented oval portholes operated by a wireless foot switch, providing the option of gloved or bare-handed work. Accommodating up to 20 plates, each of these mini-airlock portholes enables users to save time and gas when just a few samples have to be transferred.

The latest addition to Don Whitley Scientific's high performance anaerobic workstation range, the A85 facilitates fast and effective sample manipulation and equipment transfer within a secure anaerobic environment. The unique motorized gloveport system combines ultimate performance with user comfort and reliability. It also features a 30 litre airlock that accommodates up to 90 petri dishes. This ensures that only the minimum amount of oxygen is introduced into the chamber.

Created with research and clinical scientists in mind, the A85 offers optimised incubation capacity to suit small and large sample quantities. Both the working and storage areas of the chamber are temperature controlled with a uniform temperature



don whitley
scientific

gradient throughout. This allows the full internal chamber volume to be used for incubation, with a maximum capacity of up to 540 petri dishes. Every chamber is heat treated during manufacture to increase the durability of the fabricated structure and help ensure a long working life.

Don Whitley Scientific delivers maximum flexibility across the whole equipment portfolio, developing workstations that can accommodate individual needs. The A85 is available with a variety of accessories and factory-fitted options ranging from a refrigeration unit to gas-leak detector. Should electrical equipment be required inside the workstation, the A85 can be fitted with an interior electrical socket. Workstation installation, scheduled maintenance and worldwide technical support are provided through Don Whitley Scientific's extensive network of highly-skilled engineers.





don whitley
scientific

ENDS

Editor's note

For further information, please contact: Deborah Robinson, Marketing Manager,
Don Whitley Scientific, 14 Otley Road, Shipley, West Yorkshire, BD17 7SE, UK.

Tel: +44 (0)1274 595728 Fax: +44 (0)1274 531197

Email: deborah_robinson@dwscientific.co.uk

Website: www.dwscientific.co.uk

For further press information, please contact: Kat Kalinina / Nicola Aldern, Barrett
Dixon Bell, Craig Court, 25 Hale Road, Altrincham, Cheshire WA14 2EY, UK. Tel: +44
(0)161 925 4700 Fax: +44 (0)161 925 4701

Email: kat@bdb.co.uk / nicola@bdb.co.uk